Periodontal Regeneration



Treatment concept of Dr. Pierpaolo Cortellini, Florence, Italy



- > Periodontal regeneration in an aesthetic area
- > Severe generalized periodontitis
- > Severe generalized gingival inflammation and large deposits of plaque
- > Deep (10 mm) and wide 1-wall intrabony defect

1. Indication profile

marcation prome		
Region	× aesthetic region	non-aesthetic region
	single tooth gap	multiple tooth gap
Bony situation	imes bone defect present	no bone defect present
Soft tissue situation	recession	× no recession
	imes inflamed	infected
	imes thick biotype	thin biotype
	imes primary wound closure possible	primary wound closure not possible
	× intact papillae	impaired, missing papillae
	imes adequate keratinised mucosa	inadequate keratinised mucosa uneventful
Further periodontal examination	imes full mouth plaque score: 99%	
	imes full mouth bleeding score: 100%	

Background information

Patient:

Male, 34 years old, referred to the practice for periodontal therapy.

First examination: January 2002.

Chief compliants: Inflammation, recurrent abscesses, pain, migration and hypermobility of tooth 11, bad breath. The patient was concerned with function, aesthetics, and preservation of teeth.

Anamnesis: Good general health, family history of periodontitis, never smoker, never treated for periodontitis.

Periodontal examination: Severe generalized gingival inflammation associated with the presence of large deposits of plaque and calculus, migration of tooth 11, purulence associated to tooth 11 and 41. Bad breath. Full mouth plaque score: 99%. Full mouth bleeding score: 100%. Presenting with 108 sites with probing depth ≥5 mm.

Diagnosis: Chronic generalized severe periodontal disease in a patient with family history for periodontitis and presence of large deposits of plaque and calculus.

Initial treatment plan: Cause related periodontal therapy, including motivation and instructions for home care, professional supra-gingival debridement and sub-gingival root planing. Re-evaluation for potential additional therapy.

Treatment objectives: Gain periodontal health, preserve teeth, improve function and aesthetics.

Re-evaluation: 1 month after completion of cause-related therapy the patient reported the complete resolution of bad breath, resolution of inflammation and purulence, resolution of pain, lower mobility associated to tooth 11. Tooth 11 appeared also slightly repositioned with respect to baseline migration. Full mouth plaque score: 17%. Full mouth bleeding score: 10%. Presenting now with 13 sites with residual probing depth \geq 5 mm. Residual pockets were associated to teeth 16-17 and tooth 11. Radiographic examination showed the presence of a deep intrabony defect associated with tooth 11 1,2,3,4,5 .

Surgical treatment plan: Flap surgery teeth 16-17; periodontal regeneration tooth 11.

2. Aims of the therapy

> Aims of periodontal regeneration tooth 11: At re-evaluation, tooth 11 presented with residual pockets of 11 mm at the mesial side, 9 mm at the palatal side, and 7 mm at the distal side, associated with a deep and wide intrabony defect. Soft tissues were well preserved and represented by a consistent amount of thick attached gingiva. Periodontal regeneration was planned to reduce probing depth by increasing bone and attachment in order to avoid gingival recession and to reduce tooth hypermobility. Overall aims, therefore, were resolution of pockets, aesthetic preservation and function improvements^{6,7,8}.

3. Surgical procedure




Fig. 1 Baseline photograph: evidence of severe gingival inflammation, associated with plaque and calculus accumulation, and migration of tooth 11.



Fig. 2 Re-evaluation photograph: resolution of the gingival inflammation. Tooth 11 slightly repositioned (partial spontaneous resolution of the pathological migration).



Fig. 3 Tooth 11 after cause-related therapy presenting with a slight residual migration and no gingival recession. The gingiva is thick and the interdental papillae well preserved.



4 The baseline radiograph shows the presence of a deep and wide intrabony defect involving the



The flap design involves tooth 12 (distal angle) through tooth 21 (distal angle). After elevation of a full thickness buccal and lingual flap, a deep (10 mm) and wide 1-wall intrabony defect associated to tooth 11 is evident 1.2.12.13.



Fig. 10 After membrane adaptation, a split thickness incision has been performed on the buccal flap associated with a vertical releasing incision distal to tooth 12 to increase flap mobility and allow primary



Fig. 13 Primary closure maintained after 1 week, at suture removal.



Fig. 16 Re-evaluation after 6-years. The tooth 11 is completely and spontaneously realigned. The patient refers good comfort and function and is fully satisfied with aesthtics 20,21,22



Fig. 5 Pre-operatory slide showing the 11 mm mesial pocket. Local anesthesia has been delivered on the buccal and lingual area.



Fig. 8 After careful debridement and root planing, a bio-resorbable collagen barrier membrane (Geistlich Bio-Gide®) has been adapted and positioned around tooth 11 ¹⁴.



Fig. 11 Primary closure of the flap has been achieved with multilayer internal mattress sutures. The interdental space between 11 and 21 has been closed with 3 levels of sutures 12,18.



Fig. 14 Re-evaluation at 1 year. The tooth 11 has spontaneously realigned, no gingival recession has occurred, and the mobility is completely resolved. A 4 mm residual probing depth is evident, along with a 7 mm attachment level gain as compared to baseline with a good preservation of the interdental soft tissues 14,19.



Fig. 17 Radiograph taken 6-years after regeneration, showing the stability of the defect resolution.



Fig. 6 A modified papilla preservation incision has been performed between teeth 11 and 21 (wide interdental space), while a simplified papilla preservation flap has been preferred between teeth 12 and 11 (narrow interdental space) 9,10,7



Fig. 9 A deproteinized bovine bone mineral (Geistlich Bio-Oss®) is implanted to fill the intrabony defect and support the collagen barrier membrane. The use of a combined approach (barrier and filler) has been chosen in this case to support the gingival tissues in the presence of a non-supportive wide 1-wall intrabony defect ^{13,14,15,16,17}.



Fig. 12 Post-operatory radiograph, showing the intrabony defect filled by the implanted material.



Fig. 15 The 1-year radiograph shows the resolution of the intrabony component of the defect.

Literature references

- Cortellini, P., Pini-Prato, G. & Tonetti, M. (1993). Periodontal regeneration of human infrabony defects. I. Clinical Measures. Journal of Periodontology 64, 254-260.
- ² Cortellini, P., Pini-Prato, G. & Tonetti, M. (1993). Periodontal regeneration of human infrabony defects. II. Re-entry procedures and bone measures. Journal of Periodontology 64, 261268.
- Tonetti, M., Pini-Prato, G. & Cortellini, P. (1993). Periodontal regeneration of human infrabony defects. IV. Determinants of the healing response. Journal of Periodontology 64, 934-940.
- ⁴ Tonetti, M., Pini-Prato, G. & Cortellini, P. (1995). Effect of cigarette smoking on periodontal healing following GTR in infrabony defects. A preliminary retrospective study. Journal of Clinical Periodontology 22, 229-234.
- Tonetti, M., Pini-Prato, G. & Cortellini, P. (1996). Factors affecting the healing response of intrabony defects following guided tissue regeneration and access flap surgery. Journal of Clinical Periodontology 23, 548-556.
- ⁶ Cortellini, P. & Tonetti, M.S. (2000). Focus on intrabony defects: guided tissue regeneration (GTR). Periodontology 2000 22, 104-132.
- Cortellini P, Tonetti MS. (2005). Clinical performance of a regenerative strategy for intrabony defects: scientific evidence and clinical experience. J Periodontol. Mar;76(3):341-50.
- Cortellini P. Reconstructive periodontal surgery: a challenge for modern periodontology. Int Dent J 2006; 56 suppl 1:250-255.
- Cortellini, P., Carnevale, G., Sanz, M. & Tonetti, M.S. (1998). Treatment of deep and shallow intrabony defects. A multicenter randomized controlled clinical trial. Journal of Clinical Periodontology 25, 981-987.
- Cortellini, P., Pini-Prato, G. & Tonetti, M. (1996). The modified papilla preservation technique with bioresorbable barrier membranes in the treatment of intrabony defects. Case reports. International Journal of Periodontics and Restorative Dentistry 14, 8-15.
- Cortellini, P., Prato, G.P. & Tonetti, M.S. (1999). The simplified papilla preservation flap. A novel surgical approach for the management of soft tissues in regenerative procedures. International Journal of Periodontics and Restorative Dentistry 19, 589-599.
- Cortellini, P., Pini-Prato, G. & Tonetti, M. (1995). Periodontal regeneration of human infrabony defects with titanium reinforced membranes. A controlled clinical trial. Journal of Periodontology 66, 797-803.
- Cortellini, P. & Tonetti, M. (1999). Radiographic defect angle influences the outcome of GTR therapy in intrabony defects. Journal of Dental Research 78, 381 (abstract).
- Tonetti MS, Cortellini P, Lang NP, Suvan JE, Adriaens P, Dubravec D, Fonzar A, Fourmousis I, Rasperini G, Rossi R, Silvestri M, Topoll H, Wallkamm B, Zybutz M. (2004). Clinical outcomes following treatment of human intrabony defects with GTR/bone replacement material or access flap alone. A multicenter randomized controlled clinical trial. J Clin Periodontol. Sep;31(9):770-6.
- 5 Cortellini, P., Tonetti, M.S., Lang, N.P., Suvan, J.E., Zucchelli, G., Vangsted, T., Silvestri, M., Rossi, R., McClain, P., Fonzar, A., Dubravec, D. & Adriaens, P. (2001). The simplified papilla preservation flap in the regenerative treatment of deep intrabony defects: clinical outcomes and postoperative morbidity. Journal of Periodontology 72, 1701-1712.
- Linares A, Cortellini P, Lang NP, Suvan J, Tonetti MS. Guided tissue regeneration/deproteinized bovine bone mineral or papilla preservation flap alone for treatment of intrabony defects. II: radiographic predictors and outcomes. J Clin Periodontol 2006; 33: 351-358.
- Tonetti MS, Fourmousis I, Suvan J, Cortellini P, Bragger U, Lang NP; European Research Group on Periodontology (ERGOPERIO). (2004). Healing, post-operative morbidity and patient perception of outcomes following regenerative therapy of deep intrabony defects. J Clin Periodontol. Dec;31(12):1092-8.
- Cortellini, P. & Tonetti, M.S. (2001). Microsurgical approach to periodontal regeneration. Initial evaluation in a case cohort. Journal of Periodontology 72, 559-569.
- Mayfield L, Tonetti MS, Cortellini P, Lang NP. Microbial colonisation pattern predict the outcomes of surgical treatment of intrabony defects. J Clin Periodontol 2006; 33:62-68.
- 20 Cortellini, P., Pini-Prato, G. & Tonetti, M. (1994). Periodontal regeneration of human infrabony defects. V. Effect of oral hygiene on long term stability. Journal of Clinical Periodontology 21, 606-610.
- 21 Cortellini, P., Pini-Prato, G. & Tonetti, M. (1996). Long term stability of clinical attachment following guided tissue regeneration and conventional therapy. Journal of Clinical Periodontology 23, 106-111.
- 22 Cortellini P, Tonetti MS. (2004). Long-term tooth survival following regenerative treatment of intrabony defects. J Periodontol. May;75(5):672-8.

Further reading

Cortellini, P., Pini-Prato, G. & Tonetti, M. (1996). Periodontal regeneration of human intrabony defects with bioresorbable membranes. A controlled clinical trial. Journal of Periodontology 67, 217-223.

Cortellini, P., Pini-Prato, G. & Tonetti, M. (1995). The modified papilla preservation technique. A new surgical approach for interproximal regenerative procedures. Journal of Periodontology 66, 261-266.

Cortellini, P. & Tonetti, M. (2000). Evaluation of the effect of tooth vitality on regenerative outcomes in intrabony defects. Journal of Clinical Periodontology 28, 672-679.

Cortellini P, Tonetti MS (2007) A minimally invasive surgical technique (MIST) with enamel matrix derivate in the regenerative treatment of intrabony defects: a novel approach to limit morbidity. Journal of Clinical Periodontology; 34: 87-93.

Cortellini P, Tonetti MS (2007) Minimally invasive surgical technique (MIST) and enamel matrix derivative (EMD) in intrabony defects. (I) Clinical outcomes and morbidity. Journal of Clinical Periodontology; 34: 1082-1088.

Tonetti, M. S., Pini-Prato, G. P., Williams, R. C. & Cortellini, P. (1993). Periodontal regeneration of human infrabony defects. III. Diagnostic strategies to detect bone gain. Journal of Periodontology 64, 269-277.

Tonetti, M., Pini-Prato, G. & Cortellini, P. (1996). Guided tissue regeneration of deep intrabony defects in strategically important prosthetic abutments. International Journal of Periodontics and Restorative Deptistry 16, 378-387.

Tonetti, M., Cortellini, P., Suvan, J.E., Adriaens, P., Baldi, C., Dubravec, D., Fonzar, A., Fourmosis, I., Magnani, C., Muller-Campanile, V., Patroni, S., Sanz, M., Vangsted, T., Zabalegui, I., Pini Prato, G. & Lang, N.P. (1998). Generalizability of the added benefits of guided tissue regeneration in the treatment of deep intrabony defects. Evaluation in a multi-center randomized controlled clinical trial. Journal of Periodontology 69, 1183-1192.

Tonetti, M., Lang, N.P., Cortellini, P. et al. (2002). Enamel matrix proteins in the regenerative therapy of deep intrabony defects. A multicenter randomized controlled clinical trial. Journal of Clinical Periodontology;29:317-25.

Tsitoura E, Tucker R, Suvan J, Laurell L, Cortellini P, Tonetti M. Baseline radiographic defect angle of the intrabony defect as a prognostic indicator in regenerative periodontal surgery with enamel matrix derivative. J Clin Periodontol 2004; 31:643-647.

Suppliers

- > Blades: Swann-Morton #15, Swann-Morton LTD, Sheffield, England and Micro USM 6900, Sable Industries, Vista CA, USA
- > Sutures: Gore-tex CV-6 P6K23A needle RT 13 and CV-7 P7K13A needle RT 11, W.L. Gore & Ass, Flagstaff AZ, USA

Contact

Dr. Pierpaolo Cortellini, Via Carlo Botta 16, 50136 Florence, Italy telephone: +39 055 243950, fax: +39 055 2478031, e-mail: studiocortellini@cortellini.191.it

Further Indication Sheets

- > For free delivery please contact: www.geistlich.com/indicationsheets
- > If you no longer wish to collect Indication Sheets, please unsubscribe with your local distribution partner



©Geistlich Pharma AG
Business Unit Biomaterials
CH-6110 Wolhusen
phone +41 41 492 56 30
fax +41 41 492 56 39
www.geistlich.com